



Stablcal Stabilised Formazine Turbidity Standard, 1000 NTU, 500 mL

Product #: 2660649
ZAR Price: Contact Hach
Ships within 3 days

Ready to use, stabilised Formazine Primary standard solution in 500 mL bottle, 1000 NTU. For TL23, TU52, 2100Q, 2100N, 2100AN series.

Specifications

Concentration:	1000 NTU
Platform:	Solution
Quantity:	500 mL

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)