



## Calibration kit, Stablcal turbidity standards, 2100N / N IS turbidimeter, 500 mL bottles

**Product #:** 2662100  
**ZAR Price:** Contact Hach  
**Ships within 2 weeks**

Five ready to use, stabilised Formazine Primary standard solutions in bottles. Kit includes 500 mL of each: <0.1 NTU, 20 NTU, 200 NTU, 1000 NTU, 4000 NTU.

---

### Specifications

Concentration:	< 0.1, 20, 200, 1000, 4000 NTU
Quantity:	Set of 5

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)