



Secondary Turbidity Standards Kit for 2100 Portable Turbidimeters

Product #:

2464105

ZAR Price:

Contact Hach

Ships within 3 days

Kit includes 1 sealed vial each, 0-10 NTU, 10-100 NTU, 100-1000 NTU.

Gelex Secondary Standards Kit, for 2100P and 2100Q Portable Turbidimeters (includes 1 sealed vial each, 0-10 NTU, 10-100 NTU, 100-1000 NTU).

Specifications

Quantity: Set of 3

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)