



Gelex Secondary Turbidity Standards Kit, 0 - 10000 NTU

Product #:	2589200
ZAR Price:	Contact Hach
Available	

Kit includes 1 sealed vial each, Stray Light Standard, 0-2 NTU, 0-20 NTU, 0-200 NTU, 200-4000 NTU, 4000-10000 NTU. For 2100N, 2100N IS, 2100AN, 2100AN IS.

Specifications

Quantity:	Set of 6
-----------	----------

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)