



Calibration kit, Stablcal turbidity standards, 2100N / N IS turbidimeter, 100 mL bottles

Product #: 2662110
ZAR Price: Contact Hach
Ships within 3 days

Stablcal Standards Calibration Kit, for 2100N and 2100N IS Laboratory Turbidimeters. USEPA accepted for calibration purposes. Stablcal Stabilized Formazin Turbidity Standards are true Formazin dilutions developed for use in any turbidimeter. Kit includes 100 mL bottles of <0.1, 20, 200, 1000 and 4000 NTU standards.

Precise, pre-mixed Formazin solutions

With proprietary manufacturing technology, Hach prepares Stablcal Standards in precise concentrations at a wide variety of concentrations - providing accurate turbidimeter calibrations.

Based on Formazin Primary Standards

Stablcal Stabilised Formazin Standards contain the same light scattering polymer as traditional Formazin primary turbidity standards - and perform equivalently.

Enhanced stability saves time and labour

While Formazin primary standards can be prepared directly in the lab, the process is labor intensive and time consuming - and the diluted standards are not stable. Stablcal is a stabilised version of Formazin - providing repeatable results, without the time and effort required for lab-prepared standards.

Confidence in results

Hach was the first to apply Formazin as a turbidity standard in the drinking water industry, more than 40 years ago. Stablcal standards are based on decades of experience in turbidity measurement.

Specifications

Concentration:	< 0.1, 20, 200, 1000, 4000 NTU
Quantity:	Set of 5

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)