



## Stablcal Ampule Calibration Kit, Formazin, for 2100Q/QIS portable turbidimeter

<b>Product #:</b>	<b>2971205</b>
ZAR Price:	Contact Hach
Available	

Calibration Standard kit for 2100Q/QIS portable turbidimeter. Ampule kit contains Formazin standards in concentrations: 10, 20, 100, 800 NTU.

**Precise, Pre-mixed Formazin Solutions**

**Based on Formazin Primary Standards**

**Enhanced Stability Saves Time and Labor**

**Confidence in results**

---

### Specifications

Concentration:	20, 100, 800, 10 NTU (for verification)
----------------	---

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)