



## Ferric Ion Solution (100 mL)

**Product #:** 2212242  
**ZAR Price:** Contact Hach  
**Ships within 3 days**



Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Reagent solution for colorimetric Chloride determination, 0.1-25.0 mg/L  $\text{Cl}^-$ .

---

### Specifications

Method:	Mercuric Thiocyanate
Number of tests:	100
Parameter:	Chloride
Platform:	Solution
Quantity:	100 mL

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)