



Silver nitrate digital titrator cartridge, 0.2256 N, 1 piece

Product #: 1439601
ZAR Price: Contact Hach
Ships within 3 days

Ready-to-use titrant solution for the analysis of Chloride (Silver Nitrate Method), 10 - 100 mg/L Cl^- .

Approximately 40-100 tests.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifications

Method:	8207
Number of tests:	40 - 100
Parameter:	Chloride
Platform:	Digital Titrator
Range:	10 - 100 mg/L Cl^-
Titrant:	0.2256 N AgNO_3

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)