



Silver nitrate digital titrator cartridge, 1.128 N, 1 piece

Product #:

1439701

ZAR Price:

Contact Hach

Ships within 3 days

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Ready-to-use titrant solution for the analysis of Chloride (Silver Nitrate Method), 10 - 10000 mg/L Cl^- .

Approximately 40-100 tests.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifications

Method:	8207
Number of tests:	40 - 100
Parameter:	Chloride
Platform:	Digital Titrator
Range:	100 - 10000 mg/L Cl^-
Titrant:	1.128 N AgNO_3

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)