



Chloride Reagent Set, 0.1-25.0 mg/L Cl⁻

Product #:

2319800

ZAR Price:
Available

Contact Hach

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Contains mercury. Dispose of according to local, state and federal regulations.

Reagent set for the determination of Chloride, Mercuric Thiocyanate method.

Specifications

Footnote:	*Generally, the DR 5000 performs the same methods as the DR/4000. The DR 2800 performs the same methods as the DR 2700, DR/2500, and DR/2400. Some exceptions do exist.
Instrument:	All except DR 900
	DR/800 Series, PC II
Method:	Mercuric Thiocyanate
Number of tests:	100
Parameter:	Chloride
Platform:	Solution
Range:	0 - 25.0 mg/L Cl ⁻

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)