



## Chloride Reagent Set, 10-8000mg/L Cl

Product #:

2272600

ZAR Price:  
Available

Contact Hach

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

**Contains mercury. Dispose of according to local, state and federal regulations.**

For determination of chloride by mercuric nitrate Digital Titrator Method (Hach Method 8206). Reagent Set includes Diphenylcarbazone Reagent Powder Pillows, Mercuric Nitrate Digital Titrator Cartridges.

---

### Specifications

Method:	Mercuric Nitrate Titration
Number of tests:	100
Platform:	Digital Titrator
Range:	10 - 8000 mg/L Cl <sup>-</sup>

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)