



Chloride Reagent Set, 10-10000mg/L Cl

Product #:

2288000

ZAR Price:

Contact Hach

Ships within 3 days

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

For the determination of chloride by Silver Nitrate Digital Titrator Method (Hach Method 8207). Set contains Chloride 2 Indicator Powder Pillows and two Silver Nitrate Titration Cartridges.

Specifications

Method:	8207
Number of tests:	100
Platform:	Digital Titrator
Range:	10 - 10000 mg/L Cl ⁻

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)