



Cyanide Reagent Set, Powder Pillows, 0.002-0.240 mg/L CN

Product #:	2430200
ZAR Price:	Contact Hach
Available	

Powder reagents in sealed packages for the determination of Cyanide, Pyridine-Pyrazalone Method. For 10 mL sample. 100 tests.

Specifications

Method:	8027
Number of tests:	100
Parameter:	Cyanide
Platform:	Powder Pillows
Range:	0.002 - 0.240 CN ⁻
Sample size:	10 mL

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)