



Reagent powder pillows, SulfaVer 4, sulfate , 25 mL, pk/1000

Product #: 1206528

OBSOLETE ITEM
This item is no longer available.

Suggested replacements:

- 1206599

Powder reagent in sealed packages for the determination of Sulfate. For 25ml sample.

For sulfate determination by the SulfaVer 4 Turbidimetric Method. USEPA accepted for reporting wastewater analysis using Hach Method 8051. Range: up to 70 mg/L. Sample Size: 25 mL. Pack of 1000 powder pillows.

Specifications

EPA compliant:	Yes
Instrument:	DR/4000, DR/3000, DR/2010, DR/2000
Method:	8051
Number of tests:	1000
Platform:	Powder Pillows
Quantity:	1000 /pk
Range:	2 - 70 mg/L SO ₄

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)