



Sulphate Reagent Powder Pillows, 2-70 mg/L SO₄

Product #: 2106769
ZAR Price: Contact Hach
Available

Powder reagent in sealed packages for the determination of Sulphate by the SulfaVer turbidimetric method. For 10mL sample; 100/pk.

Specifications

EPA compliant:	Yes
Instrument:	All except DR/4000, DR/3000, DR/2010, DR/2000
Method:	8051
Number of tests:	100
Parameter:	Sulfate
Platform:	Powder Pillows
Quantity:	100/pk
Range:	2 - 70 mg/L SO ₄ ²⁻

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)