



Iodate-Iodide Digital Titrator Cartridge, 0.3998 N

Product #:

1496101

ZAR Price:
Available

Contact Hach

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Ready-to-use titrant solution for the analysis of Sulfite, 4 - >400mg/L SO_3 .

Approximately 40-100 tests.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifications

Method:	Iodate-Iodide
Parameter:	Sulfite
Platform:	Digital Titrator
Range:	0 - 800 mg/L SO_3^{2-}
Titrant:	0.3998 N KIO_3 -KI

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)