



## Reagent set, sulfite, SU-5 kit

**Product #:**

**2440000**

**ZAR Price:**

Contact Hach

**Ships within 3 days**

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Reagent set for determination of Sulfite by iodometric drop count titration. Replacement reagent set for SU-5 Test Kit (148002). Range: Up to 200 mg/L SO<sub>3</sub>. Set includes Sulfamic Acid Powder Pillows, Sulfite 1 Reagent Powder Pillows and Sulfite 3 Solution for approximately 100 tests.

---

### Specifications

Special Features:

Sulfite Reagent set, SU-5 Kit : Iodometric Method

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)