



## Reagent set, sulfite, buret titration

**Product #:**

**2459800**

**ZAR Price:**

Contact Hach

**Ships within 3 days**



**Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

For determination of Sulfite by iodometric buret titration. Hach Method 8071 Range: 0 to 500+ mg/L  $\text{SO}_3^{2-}$ . Reagent set includes: Dissolved Oxygen 3 Reagent Powder Pillows, Potassium Iodide-Iodate Standard Solution and Starch Indicator Solution for approximately 100 tests.

---

### Specifications

Method:	8071
Number of tests:	100
Range:	0 - > 500 mg/L $\text{SO}_3^{2-}$
Special Features:	Sulfite Reagent set, Buret Titration

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)