



## Chlorine Dioxide Reagent Set, 0.01-1.00 mg/L $\text{ClO}_2$

**Product #:** 2242300  
**ZAR Price:** Contact Hach  
**Ships within 3 days**

For the determination of Chlorine Dioxide, Chlorophenol Red Method. For 50mL sample. 100 tests.

---

### Specifications

|                  |                                            |
|------------------|--------------------------------------------|
| Instrument:      | All except DR/800 Series,<br>DR 900, PC II |
| Method:          | Chlorophenol Red                           |
| Number of tests: | 100                                        |
| Parameter:       | Chlorine Dioxide                           |
| Range:           | 0.01 - 1.00 mg/L $\text{ClO}_2$            |

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)