



Chlorine Dioxide Reagent Set, DPD Glycine Accuvac Ampules

Product #: 2771000

 **Hazardous**

OBsolete ITEM
This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

Suggested replacements:
• 2770900

If necessary, we will change your selected shipping method to accommodate these items.

For determination of Chlorine Dioxide by the DPD Glycine method, using Accuvac Ampules. USEPA accepted for reporting for drinking water analysis using Hach Method 10126. Range: 0.04 - 5.00 mg/L ClO_2 . Reagent set includes DPD Free Chlorine Accuvac Ampules and Glycine Reagent for 25 tests.

Boric Acid free!

Specifications

EPA compliant:	Yes
Instrument:	All
Method:	10126
Number of tests:	25
Parameter:	Chlorine Dioxide
Range:	0.04 - 5.00 mg/L ClO_2

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)