



DPD Free Chlorine Reagent Powder Pillows, 10 mL, pk/300

Product #:	2105503
ZAR Price:	Contact Hach
Available	

For free chlorine determination by the DPD Free Chlorine method.

USEPA accepted for reporting for drinking water analysis using Hach Method 8021.

Range: 0.02 to 2.00 mg/L Cl₂. Sample Size: 10 mL.

Pack of 300 powder pillows.

*This product has not been evaluated to test for chlorine and chloramines in medical applications in the United States.

Specifications

EPA compliant:	Yes
Instrument:	All except DR/3000, DR/2000
Method:	DPD
Number of tests:	300
Parameter:	Free Chlorine
Platform:	Powder Pillows
Range:	0.02 - 2.00 mg/L Cl ₂
Sample size:	10 mL
Storage conditions:	10 - 25 °C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)