



## Free Chlorine DPD Reagent Powder Pillows, 10 mL, pk/1000

**Product #:** 2105528  
**ZAR Price:** Contact Hach  
**Available**

### For Free Chlorine determination by the DPD method

These Free Chlorine Powder Pillow reagents are USEPA accepted for reporting in drinking water and wastewater analysis using Hach Method 8021. For best results use with a Hach colorimeter or spectrophotometer.

#### Product Details

- Pack of 1000 powder pillows
- 10 mL sample size
- Range from: 0.02 - 2.00 mg/L  $\text{Cl}_2$
- Boric Acid free!

Want to simplify reagent reordering? Ask your Hach representative about our simple and easy automatic delivery program for Hach reagents today.

---

### Specifications

EPA compliant:	Yes
Instrument:	All except DR/3000, DR/2000
Method:	8021
Number of tests:	1000
Parameter:	Free Chlorine, Chlorine Dioxide
Platform:	Powder Pillows
Range:	0.02 - 2.00 mg/L $\text{Cl}_2$
Sample size:	10 mL
Storage conditions:	10 - 25 °C

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)