



## Free Chlorine DPD Reagent Powder Pillows, 10 mL, pk/100

**Product #:** 2105569  
**ZAR Price:** Contact Hach  
**Available**

For free chlorine determination by the DPD Free Chlorine method. USEPA accepted for reporting for drinking water analysis using Hach Method 8021.

Range: 0.02 - 2.00 mg/L  $\text{Cl}_2$ . Sample Size: 10 mL.

Pack of 100 powder pillows.

Boric Acid free!

---

### Specifications

Description:	DPD Free Chlorine Reagent Powder Pillows, 10 mL
EPA compliant:	Yes
Instrument:	All except DR/3000, DR/2000
Method:	8021
Number of tests:	100
Parameter:	Free Chlorine, Chlorine Dioxide
Platform:	Powder Pillows
Range:	0.02 - 2.00 mg/L $\text{Cl}_2$
Sample size:	10 mL
Storage conditions:	10 - 25 °C

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)