



Free Chlorine Reagent Powder Pillows, DPD (pk/100)

Product #:	1407799
ZAR Price:	Contact Hach
Available	

For determination of free chlorine by the DPD method. For use with Hach's Chlorine Visual Test Kits.

Sample Size: 5 mL. Range: up to 3.5 mg/L Cl_2 .

Pack of 100 powder pillows.

Boric Acid free!

Specifications

Method:	DPD
Number of tests:	100
Parameter:	Chlorine, Free
Platform:	Powder Pillows
Sample size:	5 mL
Storage conditions:	10 - 25 °C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)