



DPD Free Chlorine Reagent Powder Pillows, 5 mL, pk/50

Product #: 2197846
ZAR Price: Contact Hach
Ships within 3 days

For determination of free chlorine with Hach's Chlorine Visual Test Kits- 2060300. Sample Size: 5 mL. Range: up to 2.5 mg/L Cl_2 . Pack of 50 Powder Pillows.

Boric Acid free!

Specifications

Method:	DPD
Number of tests:	50
Parameter:	Chlorine, Free
Platform:	Powder Pillows
Sample size:	5 mL
Storage conditions:	10 - 25 °C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)