



DPD Free chlorine, Swiftest dispenser reagent (refill)

Product #:	2105560
ZAR Price:	Contact Hach
Available	

For free chlorine determination by the DPD Free Chlorine method. For use with Hach Swiftest DPD dispenser.

USEPA accepted for reporting for drinking water analysis using Hach Method 8021. Total chlorine and potassium iodide refill vials also available.

Range: 0.02 - 2.00 mg/L Cl_2 . Sample size: 10 mL. Approximately 250 tests per vial.

*This product has not been evaluated to test for chlorine and chloramines in medical applications in the United States.

Specifications

EPA compliant:	Yes
Method:	DPD
Number of tests:	250
Parameter:	Chlorine, Free
Platform:	Swiftest
Range:	0.02 - 2.00 mg/L Cl_2

Required Accessories

- Swiftest DPD Free chlorine reagent dispenser and reagent vial (Item 2802300)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)