



## Swiftest DPD Free chlorine reagent dispenser and reagent vial

**Product #:** 2802300  
**ZAR Price:** Contact Hach  
**Ships within 3 days**

For the determination of Free Chlorine by the DPD method, 0.02-2.00 mg/L Cl<sub>2</sub>.

---

### Specifications

Instrument:	All except DR/3000, DR/2000
Method:	DPD
Number of tests:	Approximately 250
Parameter:	Chlorine, Free
Platform:	Swiftest
Range:	0.02 - 2.00 mg/L Cl <sub>2</sub>

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)