



Free Chlorine Reagent, Accuvac Ampules, DPD, pk/25

Product #:	2502025
ZAR Price:	Contact Hach
Ships within 3 days	

Ampules contain the precise amount of reagent for a single test and can be used as a measurement cuvette. For the determination of Free Chlorine, 0.02-2.00 mg/L Cl_2 .

Specifications

EPA compliant:	Yes
Method:	DPD
Method Number 2:	10126
Number of tests:	25
Parameter:	Free Chlorine, Chlorine Dioxide
Platform:	Accuvac
Quantity:	25 /pk
Range:	0.02 - 2.00 mg/L Cl_2

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)