



## DPD Total Chlorine Reagent Powder Pillows, 5 mL, pk/100

**Product #:** 1407699  
**ZAR Price:** Contact Hach  
**Ships within 3 days**

For determination of total chlorine by the DPD method. For use with Hach's Chlorine Visual Test Kits. Sample Size: 5 mL. Range: up to 3.5 mg/L  $\text{Cl}_2$ . Pack of 100 powder pillows.

Boric Acid free!

---

### Specifications

Method:	DPD
Number of tests:	100
Platform:	Powder Pillows
Sample size:	5 mL
Storage conditions:	10 - 25 °C

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)