



DPD Total chlorine, Swiftest dispenser reagent (refill)

Product #: 2105660
ZAR Price: Contact Hach
Ships within 3 days

DPD Reagent dispenser replacement vial for the determination of Total Chlorine, 0.02-2.00 mg/L Cl₂.

Approximately 250 tests

Specifications

EPA compliant:	Yes
Instrument:	All except DR/3000, DR/2000
Method:	DPD
Number of tests:	250
Parameter:	Chlorine, Total
Platform:	Swiftest
Range:	0.02 - 2.00 mg/L Cl ₂
Range 2:	0.02 - 4.50 mg/L Cl ₂

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)