



## Swiftest DPD Total chlorine reagent dispenser and reagent vial

|                   |                |
|-------------------|----------------|
| <b>Product #:</b> | <b>2802400</b> |
| ZAR Price:        | Contact Hach   |
| Available         |                |

For the determination of Total Chlorine, 0.02-2.00 mg/L Cl<sub>2</sub> and 0.02-4.50 mg/L Cl<sub>2</sub>.

---

### Specifications

|                  |                                  |
|------------------|----------------------------------|
| Instrument:      | All except DR/3000, DR/2000      |
| Method:          | 8167                             |
| Number of tests: | Approximately 250                |
| Parameter:       | Chlorine, Total                  |
| Platform:        | Swiftest                         |
| Range:           | 0.02 - 2.00 mg/L Cl <sub>2</sub> |
| Range 2:         | 0.02 - 4.50 mg/L Cl <sub>2</sub> |

---

### Required Accessories

- DPD Total chlorine, Swiftest dispenser reagent (refill) (Item 2105660)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)