



**Be Right™**



## Total Chlorine Chemkey Reagents (pack of 300)

|                     |                |
|---------------------|----------------|
| <b>Product #:</b>   | <b>8499400</b> |
| ZAR Price:          | Contact Hach   |
| Ships within 1 week |                |

Total Chlorine Chemkeys, range 0.04 - 10 mg/L  $\text{Cl}_2$ . For SL1000 Portable Parallel Analyser (PPA). 12 boxes of 25 (300 tests total quantity).

---

### Specifications

|             |                              |
|-------------|------------------------------|
| Instrument: | SL1000                       |
| Method:     | 10260                        |
| Parameter:  | Chlorine, Total              |
| Platform:   | Chemkeys                     |
| Quantity:   | pk/300                       |
| Range:      | 0.04 - 10 mg/L $\text{Cl}_2$ |