



Total Chlorine Reagent, Accuvac Ampules, DPD, pk/25

Product #: 2503025
ZAR Price: Contact Hach
Ships within 3 days

Ampules contain the precise amount of reagent for a single test and can be used as a measurement cuvette. For the determination of Total Chlorine, 0.02-2.00 mg/L Cl_2 , Bromine, 0.05-4.50 mg/L Br_2 and Iodine, 0.07-7.00 mg/L I_2 .

Specifications

EPA compliant:	Yes
Instrument:	All
Method:	8167 Chlorine
Method Number 2:	8016 Bromine
Method Number 3:	8031 Jod
Number of tests:	25
Parameter:	Total Chlorine, Bromine, Iodine
Platform:	Accuvac
Quantity:	25 /pk
Range:	0.02 - 2.00 mg/L Cl_2
Range 2:	0.05 - 4.50 mg/L Br_2
Range 3:	0.07 - 7.00 mg/L I_2

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)