



## Potassium iodide powder pillows for test kits and lab methods, 100 pcs

**Product #:** 107799  
**ZAR Price:** Contact Hach  
**Ships within 3 days**

Powder reagent in sealed packages for the determination of total chlorine with lab titration and chlorine reagent set (drop count).

Reagent for several iodometric titration methods. Pack of 100 powder pillows.

---

### Specifications

Method:	10025
Method Number 2:	10026
Method Number 3:	8161
Method Number 4:	8209
Number of tests:	100
Platform:	Powder Pillows
Quantity:	100 /pk
Range:	Various

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)