



## Potassium Iodide Swiftest Bulk Dispenser, with Potassium Iodide Refill Vial

**Product #:** 2834100  
**ZAR Price:** Contact Hach  
**Ships within** 3-5 weeks

For total chlorine determination with the AT1000 Chlorine Amperometric Titrator. Includes Swiftest dispenser and vial containing enough reagent for approximately 300 tests.

### Fast, easy reagent addition

Ideal for high volume total chlorine testing using Hach's high quality potassium iodide powder.

### Just squeeze to add powder

No need to open reagent containers, less packaging waste

### Color coded by parameter

The handle of the Swiftest dispenser is color-coded for easy parameter recognition.

---

## Specifications

Method:	Amperometric Titration
Platform:	Swiftest

---

## Required Accessories

- Potassium Iodide Swiftest Dispenser Refill (Item 107760)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)