



Sodium Thiosulfate Standard Solution, Stabilised, 0.0246 N, 100 mL MDB

Product #: 2409232
ZAR Price: Contact Hach
Ships within 3 days

Stabilised sodium thiosulfate titrant for chlorine determination by iodometric titration. 100 mL Marked Dropping Bottle.

Specifications

Method:	Iodometric Titration
Number of tests:	Varies
Platform:	Solution

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)