



Freechlor F Reagent Solution, 50 mL

Product #:	2964926
ZAR Price:	Contact Hach
Ships within 3 days	

Used to determine free chlorine in the presence of manganese, chloramines and other interfering oxidants. Hach Method 10241 Requires Monochlor F Reagent, PN 2802299. 50 mL bottle. Approximately 100 tests.

Specifications

Method:	Freechlor F
Platform:	Solution

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)