



Buffer, Phosphate Type, pH 7, for Amperometric Chlorine, 100 mL MDB

Product #: 2155332
ZAR Price: Contact Hach
Ships within 3 days

For determination of chlorine by amperometric titration. APHA formulation 100 mL bottle

Specifications

Instrument:	AutoCat™ 9000
Method:	8334
Number of tests:	100
Platform:	Solution
Range:	0.5 mg/L and above Cl ₂
Range 2:	0.0051 - 5.000 mg/L Cl ₂

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)