



Phenylarsine Oxide (PAO) Standard Solution, 0.00564 N, 1 L

Product #:	199953
ZAR Price:	Contact Hach
Ships within 3 days	

For chlorine determination by amperometric titration methods. Note: also available in Digital Titrator cartridges. 0.00564 N, 1 liter bottle

Specifications

Instrument:	AutoCat™ 9000
Method:	8334
Method Number 2:	10025
Method Number 3:	8168
Platform:	Solution
Range:	0.5 mg/ L and above Cl_2
Range 2:	0.0051 - 5.000 mg/L Cl_2
Range 3:	0.0012 - 5.000 mg/L Cl_2

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)