



Acetate Buffer Solution, pH 4.0, 100 mL

Product #: 1490932
ZAR Price: Contact Hach
Ships within 3 days

Acetate Buffer Solution, pH 4.0. Total and/or Free Chlorine. Method: Amperometric. pk/100 mL.

Specifications

Instrument:	AutoCat™ 9000
Method:	Amperometric Titration
Method Number 2:	8168
Method Number 3:	10026
Method Number 4:	8209
Number of tests:	100
Platform:	Solution
Range:	0.0051 - 5.000 mg/L Cl ₂
Range 2:	0.0012 - 5.000 mg/L Cl ₂
Range 3:	6 - 1000 µg/L Cl ₂
Range 4:	15 - 1000 µg/L Cl ₂

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)