



## Fluoride Ionic Strength Adjustor (ISA) Powder Pillows, pk/100

**Product #:** 258999  
**ZAR Price:** Contact Hach  
**Ships within 3 days**

Ionic Strength Adjustor (ISA) for Fluoride determinations by Ion Selective Electrode (ISE) method. Powder form prevent from diluting samples. Ideal for low Fluoride level readings. Eliminate most interferences. Pack of 100 powder pillows.

---

### Specifications

Concentration:	ISA
EPA compliant:	Yes
Package Type:	Powder Pillows
Parameter:	Fluoride
Quantity:	100 /pk
Range:	0.01 mg/L to saturated

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)