



Fluoride Ionic Strength Adjustor (ISA), 454 g

Product #:	258901
ZAR Price:	Contact Hach
Available	

Fluoride Total Ionic Strength Adjustor Buffer (TISAB) for use in fluoride determinations by electrode methods. Eliminates interferences from aluminum and iron. Powder form does not dilute samples. Ideal for determining very low fluoride levels. 454 g jar.

Specifications

EPA compliant:	Yes
Platform:	Jar
Range:	0.01 mg/L to saturated

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)