



Fluoride reagent solution, 0.02-2.00 mg/L F (500 mL)

Product #:

44449

ZAR Price:
Available

Contact Hach

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Specifications

EPA compliant:	Yes
Instrument:	All except DR/820
Method:	SPADNS
Number of tests:	125
Parameter:	Fluoride
Platform:	Solution
Quantity:	500
Range:	0.02 - 2.00 mg/L F ⁻

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)