



Fluoride Ionic Strength Adjustor (ISA), 3.78 L

Product #:	2829017
ZAR Price:	Contact Hach
Available	

Fluoride Total Ionic Strength Adjustor Buffer Solution (TISAB) for fluoride determinations by electrode methods. Eliminates interferences from aluminum and iron.

Specifications

EPA compliant:	Yes
Number of tests:	750
Platform:	Solution
Range:	0.01 mg/L to saturated

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)