



SPADNS 2 (Arsenic-free) Fluoride Reagent Solution, 3.78 L

Product #:

2947517

ZAR Price:
Available

Contact Hach

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

For determination of Fluoride by Hach SPADNS 2 Method 10225. The SPADNS2 reagent formulation removes interference from chlorine without arsenic. Range: 0.02 - 2.00 mg/L F⁻. EPA compliant. 3.78 liter bottle

Identical performance to traditional Hach SPADNS

Use the same Hach instruments, test procedures, and calibration curves

Specifications

EPA compliant:	Yes
Instrument:	All except DR 900
Method:	SPADNS 2
Number of tests:	945
Platform:	Solution
Range:	0.02 - 2.00 mg/L F ⁻

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)