



SPADNS 2 (Arsenic-free) Fluoride Reagent Solution, 500 mL

Product #:

2947549

ZAR Price:

Contact Hach

No availability shown



Hazardous

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

For determination of Fluoride by Hach SPADNS 2 Method 10225. The SPADNS2 reagent formulation removes interference from chlorine without arsenic. Range: 0.02 - 2.00 mg/L F⁻. EPA compliant. 500 mL bottle

Identical performance to traditional Hach SPADNS

Use the same Hach instruments, test procedures, and calibration curves

Specifications

EPA compliant:	Yes
Method:	10225
Number of tests:	125
Platform:	Solution
Range:	0.02 - 2.00 mg/L F ⁻

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)