



SPADNS 2 (Arsenic-free) Fluoride Reagent Solution, 1 L

Product #:

2947553

ZAR Price:

Contact Hach

No availability shown

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

For determination of Fluoride by Hach SPADNS 2 Method 10225. The SPADNS2 reagent formulation removes interference from chlorine without arsenic. Range: 0.02 - 2.00 mg/L F⁻. EPA compliant. 1 Liter bottle.

Identical performance to traditional Hach SPADNS

Use the same Hach instruments, test procedures, and calibration curves

Specifications

| | |
|------------------|---------------------------------|
| EPA compliant: | Yes |
| Method: | 10225 |
| Number of tests: | 250 |
| Parameter: | Fluoride |
| Platform: | Solution |
| Range: | 0.02 - 2.00 mg/L F ⁻ |

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)