



SPADNS 2 Fluoride Reagent, arsenic-free, ACCUVAC, 0.02-2.00 mg/L F

Product #:
ZAR Price:
Ships within 3 days

2527025
Contact Hach



Hazardous
Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

ACCUVAC ampules contain the precise amount of reagent for a single test and can be used as a measurement cuvette. For the determination of Fluoride by the SPADNS 2 method. pk/25

Specifications

EPA compliant:	Yes
Method:	SPADNS 2
Method Number 1 Chemistry:	10225
Number of tests:	25
Parameter:	Fluoride
Platform:	Accuvac
Quantity:	25 /pk
Range:	0.02 - 2.00 mg/L F ⁻

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)