



## Ozone Accuvac, LR, 0.01-0.25 mg/L O<sub>3</sub>, pk/25

**Product #:** 2516025  
**ZAR Price:** Contact Hach  
**Ships within 3 days**

Accuvac ampules contain the precise amount of reagent for a single test and can be used as a measurement cuvette. For the determination of ozone by the Indigo method. pk/25

---

### Specifications

Description:	Ozone Accuvac, LR
Instrument:	All except DR/820
Method:	Indigo
Number of tests:	25
Parameter:	Ozone
Platform:	Accuvac
Range:	0.01 - 0.25 mg/L O <sub>3</sub>

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)