



Ozone Accuvac, 0.01-1.5 mg/L O₃

Product #: 2518025
ZAR Price: Contact Hach
Ships within 3 days

Accuvac ampules contain the precise amount of reagent for a single test and can be used as a measurement cuvette. For the determination of ozone by the Indigo method. pk/25

Specifications

Instrument:	All except DR/820, PC II
Method:	Indigo
Number of tests:	25
Parameter:	Ozone
Platform:	Accuvac
Range:	0.01 - 1.50 mg/L O ₃

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)